

Technical data

<i>Technical data</i>	<i>MINORE® 0</i>	<i>MINORE® I</i>	<i>MINORE® II</i>
<i>Model and part no. (air cooled version)</i>	MINORE® 0-A-RB400 1-100856	MINORE® I-A-RB400 1-100857	MINORE® II-A-RB400 1-100858
<i>Model and part no. (water cooled version)</i>	MINORE® 0-W-RB400 2-100271	MINORE® I-W-RB400 2-100272	MINORE® II-W-RB400 2-100273
<i>Temperature range</i>	16-22° C	16-22° C	16-22° C
<i>Cooling power @ 20°C</i>	300 Watt	350 Watt	500 Watt
<i>Pump capacity max.</i>	28 l/min.	28 l/min.	28 l/min.
<i>Feed pressure max.</i>	0.4 bar	0.4 bar	0.4 bar
<i>Connections</i>	Quick lock with 9 mm hose sleeves	Quick lock with 9 mm hose sleeves	Quick lock with 9 mm hose sleeves
<i>Dimensions W x D x H</i>	360 x 470 x 590 mm	360 x 470 x 590 mm	360 x 470 x 590 mm
<i>Current</i>	230 V/50 Hz/1 PH	230 V/50 Hz/1 PH	230 V/50 Hz/1 PH
<i>Power consumption max.</i>	220 Watt	270 Watt	340 Watt
<i>Weight</i>	31 kg	37 kg	38 kg
<i>Refrigerant</i>	R290	R290	R290